

# PROVEN PRODUCTS

## DESIGNED BY PET MOBILITY EXPERTS

Built by practitioners working on the front lines

## PROVEN EFFECTIVE

Tested in the field to ensure they'll work for you

## MEDICAL-GRADE QUALITY

Made with pet-proof materials you can trust

## OFF-THE-SHELF AFFORDABILITY, ADJUSTABILITY & SIMPLICITY

Get it faster and know they'll be able to use it at home



Walkin' Pets  
by handicappedpets.com

## PRO SERIES MOBILITY KIT



### Pro Series Mobility Kit

The Starter Kit That Gets You Everything You Need To Get Moving.

The PRO Series Mobility kit allows you to properly fit patients, no matter the condition or the size of the pet. Use this starter kit for in-house product training to educate the staff on how to use the products with patients. You can incorporate the PRO kit into rehab sessions or use to demonstrate to your customers how to use the product at home. This specialized kit makes it convenient for on the spot fitting and getting the pets back in motion.

**NOTE:** Replacement products can be purchased for 50% off. This kit also qualifies your practice for the Pro Series Referral Program.

#### This kit includes:

- 1 vertebraVe back support
- 1 Hip-EEZ base system (brace without attachments)
- 1 Front No-Knuckling Training Sock
- 1 Rear No-Knuckling Training Sock
- 1 Walkin' Fit Adjustable Splint

Walkin' Pets  
by handicappedpets.com

105 Route 101A, Suite 18 | Amherst, NH 03031

[www.walkinpets.com](http://www.walkinpets.com)

[VetResources@WalkinPets.com](mailto:VetResources@WalkinPets.com) | 866-578-2926

Walkin' Pets  
by handicappedpets.com

YOUR PROVEN PARTNER IN PET MOBILITY  
**TOGETHER WE MOVE PETS**



[www.vets.handicappedpets.com](http://www.vets.handicappedpets.com)



## Walkin' Wheels Wheelchair

The fully adjustable Walkin' Wheels wheelchair provides a simple solution for pets experiencing mobility challenges. A life-changing tool, enabling pets to exercise, play, and regain their independence. Walkin' Wheels is ideal for both long term use to keep dogs active or short-term rehabilitative use.

### Four sizes to fit ANY size pet:

Mini 2-10 lbs. | Small 11-25 lbs. | Medium 26-69 lbs. | Large 70-180 lbs

### Ideal for the following conditions:

- Hind end weakness
- Degenerative Myelopathy (DM)
- Hip Dysplasia
- Arthritis
- Paralysis
- Intervertebral Disc Disease (IVDD)
- Neurological Conditions
- Surgical or Injury recovery

### Clinical Advantages:

- Great fit with simple measurements
- Accommodates a range of pets: Fully adjustable in height, length & width
- Wheelchair adapts to patient's changing needs, easily converts into full support wheelchair
- Incorporate into physical therapy & strength training exercises



## Walkin' Fit Adjustable Splint

Brace Yourself. This Splint Does It All.

Designed to provide a comfortable, custom fit while providing the stability and support to the Tarsal or Carpal Joints. Featuring an adjustable angle to accommodate a pet's front or rear leg.

The splint's width and angle can be adjusted as you progress through a treatment plan to maintain a proper fit for the patient. Eliminates the need to stock multiple splints or the delay of custom ordering braces for different sized dogs and conditions.

### Clinical Advantages:

- Supports the Tarsal or Carpal Joints
- Patient specific anatomic considerations – "posty" or "low" conformation
- Can re-establish ROM slowly
- CCT Mechanism injury assistance
- Adjustable width accommodates bandaged limbs
- Useful for Tibiotarsal instability/Collateral ligament injuries – either medically managed or in recovery/early mobilization phase post op



## No-Knuckling Training Socks

For Front or Rear Paw Knuckling.



The No-Knuckling Training Sock is a temporary training tool, for short-term, multiple use to help correct gait and improve paw placement. Designed to enhance proprioception in patients with neurological deficits. The sensory stimulation acts as a reminder for pets to lift their limb higher during gait training.

### Clinical Advantages:

- Utilize as a gait retraining tool
- Lightweight to target patients with limb weakness
- Designed to be used by rehab professionals as part of a treatment plan
- Instruct pet owners for use during at-home training sessions
- Use During: Enhanced Cavaletti Poles, Enhanced Figure 8's or Cone Weaving, & Treadmill Patterning



## Walkin' Hip-EEZ

The First 5-in-1 Hip Support System.

The Hip-EEZ Support System is designed to provide support to the hip joint of pets. Combining the use of compression and support to alleviate hip pain. Use to treat the following conditions: Arthritis, Hip Dysplasia, Conformational Deformities, Hip Trauma Recovery, and Post-Surgical Recovery.

### Customize your patient's treatment plan by combining the Hip-EEZ with one of four attachments:

- Walkin' Hip-EEZ Bridge
  - A non-surgical Hip Subluxation solution
- Walkin' Hip-EEZ Donut
  - Decubital Ulcer prevention and treatment
- Walkin' Hip-EEZ Cross-Assist
  - Corrective training tool to prevent hind leg crossing
- Walkin' Hip-EEZ Amputee Sleeve
  - The first hip support system designed to accommodate and support the hips of rear leg amputees



## Walkin' vertebraVe

Back Support for EVERY Dog.

Give your patients the back support they need. This innovative back system for pets combines the support and comfort in a brace specifically designed for pets suffering from IVDD, arthritis and back injuries.

The dual layered support system is strategically created with a layer of medical grade memory foam that conforms to the pet's body and a second layer of 3 additional aluminum support rods placed directly over the spine for a custom fit.

### Clinical Advantages:

- Adjustable and Easy to Use
- May be used for post-surgical recovery and rehab
- Remove support rods to lessen support to increase range of motion as spine heals
- Durable and Lightweight
- Breathable material
- D-Ring for leash